

CSSD & OR
Logistical solutions

MedStor[®]
The heart of your storage solutions



Healthcare storage, logistics
and product handling





Product innovation for CSSD & OR

We specialise in the latest healthcare storage, logistics and product handling.

Instrument baskets

Our Sterisystem instruments baskets are manufactured from high polished AISI 304 graded stainless steel complete with a perfectly finished top wire to reduce the chances of damaged wraps.

2mm holes guarantee a permeability of more than 50% which benefits the washing process. The lid designed for the instrument baskets contains minimal number of components making them ideal for sterile environments.

- 2mm hole design
- Basket lid
- Colour coding system
- Silicone corner protectors



Tracking system

The Steri-ID standalone or built-in RFID tags allow for a cost-efficient, time-saving and environmentally friendly replacement of traditional instrument bar code stickers. Tags can be read with the same hand held reader used for traditional barcodes.

The Steri-ID tag is embedded in a protective silicone housing which is resistant to autoclave sterilisation conditions and also to the most common chemicals used in hospitals. The silicone housings are available in a variety of colours.

- Alternative to barcode stickers
- Protective outer coating
- Selection of colour trackers



Wire baskets & shelves

Wire baskets are high polished stainless steel, available in both Euro and DIN sizes. Extra safe and hygienic baskets are stackable and nest together making them more user-friendly. The Sterisystem high polished wire trays are the perfect solution if you need to sterilise, transport and store larger sterilisation packs.

- Stainless steel
- Available in 2 sizes
- Hygienic transportation



> Height adjustable from 765 - 1265mm

Packing tables

Ergonomically designed height adjustable table is designed to be adjusted to the task at hand. Modular accessories in addition to the tables flexibility allow you to make optimal use of the work surface.

- Ideal for users of different heights
- Easily adjustable
- Available in various sizes and finishes

Available in various sizes and finishes (HEPL, Stainless steel or HP laminate).

Specialist trolleys and mobile carts

Ideal for use in sterile environments such as Operating rooms.

Stericart

Stericart are manufactured from high polished AISI 304 graded stainless steel. They are used to move goods between the CSSD and the Operating room.

The solid closed base plate ensures instrument sets are protected from upward contamination.

- Fitted with top quality washing tunnel resistant castors
- Low rolling resistance and high manoeuvrability
- True ergonomic asset to your CSSD and OR



UN3291 aluminium transport trolleys

The closed aluminium transport trolley allow you to safely and hygienically transport sterile, non-sterile goods or medication.

Available in different widths for Euronorm and DIN baskets and shelves. The runners are mounted on easily removable ladders so the trolley can be cleaned to perfection. Steri-ID improves the tracking of containers through the logistical flow.

- Easily cleaned
- High volume transportation
- Use the Steri-ID system
- UN3291 approved

UN3291 stainless steel transport trolleys

Fitted with a lock with breakable seal and bumpers at the top and on the bottom, this transport trolley has been designed to transport sterile sets between external sterilisation unit around the hospital.

- Safely transport sterile or contaminated sets
- Washing tunnel resistant stainless steel
- Manufactured to Euronorm or DIN standards
- UN3291 approved



> Our most hygienic storage system

UFLEX storage system

This unique storage system is manufactured from anodised aluminium. The construction is completely devoid of welds so you are guaranteed a hygienic unit.

To fit out and organise your storage area the Uflex 3 storage system has been designed to accommodate a combination of Euronorm and DIN baskets and shelves.

Racks of equal depth can be added in infinity making the best possible use of the space available.

- No welds for hygienic cleaning
- Fits a combination of basket and shelves
- Available flat pack for easy installation and delivery

MedStor[®]

The heart of your storage solutions



All pictures ©Belintra NV - www.belintra.com

Mailbox
Stamford Products Limited

UK:
Bayley Street, Stalybridge, Cheshire, SK15 1QQ.
Tel: 0161 330 5577 Fax: 0161 330 5576
Website: www.medstor.co.uk email: health@medstor.co.uk

Ireland:
Unit 120, Grange Hill, Baldoyle Industrial Estate, Dublin 13.
Tel: 01 839 3596 Fax: 01 839 3468
Website: www.medstor.ie email: sales@medstor.ie