

# Procedure Trolleys & Ward Equipment

**MedStor**<sup>®</sup>  
The heart of your storage solutions



Healthcare storage, logistics  
and product handling



Dental  
Medical  
Laboratory  
Pharmaceutical

# Total space utilisation

Hygienic

Lightweight and extremely strong

Best in class

Simple to use and reconfigure

## Built around our storage systems

Conforming to International standards, both trays and dividers are injection moulded in tough clear polycarbonate, ideally suited to the Healthcare sector.

Extremely strong, robust and light weight the Medstor trays and dividers are above all designed to protect their contents.

High temperature washing machine compliant material is used to help provide hygienic surfaces within your working environment.

Simple slot in reconfigurable dividers allow trays to be both length and width ways compartmentalised in seconds - providing a comfortable and safe fit to all supplies.

Trays are also available in a grey ABS material upon request.

## Faster product Identification

Label holders are essential to increasing visibility of storage tray contents and stock control level management. All Medstor trays and divider sets are provided with a selection of label holders and coloured inserts.

## Smartkanban with RFID

The materials management system via IT technologies and bar coding utilises the Medstor labelling system and storage products to help stock control via the tried and tested 'Dual bin' system.

Please contact sales for additional information.



# Contents

# Pages

Working with you	4
Associations and accreditation	5
Ward use trolleys	6 - 13
Stainless steel dressing and instrument trolleys	14 - 15
I.V. poles and drip stands	16
Recharge trolley	17
Medical records and case note trolleys	18
General purpose	19
Linen and waste	20 - 22
Cylinder holders	23

## Introduction

Medstor is focused on delivering the best possible storage, logistics and product handling solutions to the healthcare sector.

As part of Stamford Products Ltd we take pride in our commitment to both service and quality, having already established a reputation for being market leaders in both the UK and Ireland.

**We look forward to hearing about your project requirements.**

## Product guide

For our full systems catalogue featuring all Medstor products please contact our sales office.

## Contact details

### UK:

Bayley Street, Stalybridge, Cheshire, SK15 1QQ.

Tel: 0161 330 5577 Fax: 0161 330 5576

Website: [www.medstor.co.uk](http://www.medstor.co.uk) email: [health@medstor.co.uk](mailto:health@medstor.co.uk)

### Ireland:

Unit 120, Grange Hill, Baldoyle Industrial Estate, Dublin 13.

Tel: 01 839 3596 Fax: 01 839 3468

Website: [www.medstor.ie](http://www.medstor.ie) email: [sales@medstor.ie](mailto:sales@medstor.ie)

**peto**.co.uk

Medstor products are also available via the PETO online ordering system.

Please add VO1ZZ to the start of the product code.

E.G. [VO1ZZMH1103-0000](#)





## Working with you

Analysis and consultation  
Survey and planning  
Design and costings  
Installation and implementation  
Training and education

### Every step of the way:

Our national healthcare representatives begin by assessing your current situation in terms of stock holdings, space utilisation and working practices. After consultation we will develop a tailored, individual plan based on our previous healthcare and customer experiences.

We offer a free site survey and planning service from our very own in-house team. Personalised layout plans and 3D visuals can be produced of your project including a fully itemised Medstor product listed quotation.

From a single room to a complete hospital, we work with you from start to finish for the full installation. Your project co-ordinator will liaise with you to make sure that products are delivered to your exact requirements.

Our CSCS qualified installation teams will ensure that your systems are fitted, cleaned and ready for use.

## Making specification easy

All trolleys listed feature the most efficient use of tray and divider combinations.

# Associations & accreditations

- ISO 3394: 2012 Complete, filled transport packages and unit loads
- ISO 9001: 2008 Quality Management System
- BS8555: 2003 Phase 1,2 and 3 Environmental Management System
- FSC: Forest Stewardship Council
- CHAS: Contractors Health and Safety Assessment Scheme
- BPF British Plastics Federation



# Environmental impact

Healthcare - Stamford Products Limited are committed to minimising the environmental impact of our operations in a way, which protects and enhances the environment. The company will comply with all the relevant legislation and adhere to the principles of the prevention of pollution during design and operational performance as part of our business plans.

The EMS arrangements to reduce our overall environmental impact include the following :-

- Awareness and use of appropriate technology
- Conservation of energy and reduce consumption of materials in our operations
- Reuse rather than dispose
- Encourage and promote recycling and the use of recycled and recyclable materials
- Develop products to make efficient use of materials from sustainable sources
- Establish targets and objectives to measure our improvements.



# Ward use trolleys

Sleek and stylish

Multi-functional

Solid and robust

10 year guarantee

**COMPACT TRAY 400 x 300mm CONTAINS:**

1 lengthway and widthway divider.

2 label holders and 10 coloured inserts in 5 colours.

**EURO TRAY 600 x 400mm CONTAINS:**

2 lengthway and widthway dividers.

2 label holders and 10 coloured inserts in 5 colours.

**FRAME COLOUR OPTIONS:** Lilac available from stock



Lilac



Pale green



White



Grey



Red

**CHOOSE FROM:**

Compact 400 x 300mm or Euro 600 x 400mm

**4 HEIGHT OPTIONS:**

Low level 850mm

Standard 950mm

Standard plus handles 1050mm

Tall 1450mm



KEY FEATURES:



**MEDSTOR WARD TROLLEYS CAN BE CONFIGURED AND ACCESSORISED TO SUIT ALL APPLICATIONS.**

- A** One piece mild steel welded frame with a selection of colour options. No site assembly required.
- B** Fitted with 4 anti-static 100mm castors - Non marking, buffered wheels with locking brakes.
- C** Medstor modular system trays supplied with length way and width way dividers -  
Compact tray 400 x 300mm 1 off each.  
Euro tray 600 x 400mm 2 off each.  
Clear crystal, shatterproof Medstor polycarbonate trays, available in three depths 50mm, 100mm and 200mm (max load per tray 15kgs).
- D** All units are supplied with a tough recessed seamless acrylic capped ABS top. Optional stainless steel tops are also available.
- E** Injection moulded liner sides include arresters to prevent accidental tray withdrawal. Sides and top can be quickly removed for easy cleaning.
- F** Rails on both side of the trolley to hold accessories.
- G** Lockable doors and back panels available as optional extra.

ACCESSORIES:



Ref	Stainless steel	Painted colour	Description
1	MM4479CSS	MM4479CP	90-130mm Adjustable glove box Dim. 90 x 110 x 110mm.
2	MM4479DSS	MM4479DP	Sharps bin holder - Up to 7 litre capacity Dim. 125 x 170 x 100mm.
3	MM4479ESS	MM4479EP	Gel holder Dim. 90 x 110 x 110mm.
4	MM4479FSS	MM4479FP	File holder - Accepts A4 soft and hard back files. Dim. 250 x 300 x 100mm.
5	MM4479GSS	MM4479GP	Drop down shelf - Hinged work surface extension Fits to right hand side of the trolley. 380 x 290mm.



## Low level - 850mm high

Fits under workbenches  
Ideal for seated procedures  
10 year guarantee

### OPTIONAL ACCESSORIES SEE PAGE 7

#### FOR PETO ONLINE ORDERING SYSTEM:

Please add VO1ZZ to the start of the product code.  
E.G. VO1ZZMH1103-0000

#### COMPACT TRAY 400 x 300mm CONTAINS:

1 lengthway and widthway divider.  
2 label holders and 10 coloured inserts in 5 colours.

#### EURO TRAY 600 x 400mm CONTAINS:

2 lengthway and widthway dividers.  
2 label holders and 10 coloured inserts in 5 colours.

#### FRAME COLOUR OPTIONS: Lilac available from stock



Lilac    Pale green    White    Grey    Red



Designed to store under workbench surfaces

### LOW LEVEL COMPACT SINGLE - 468 x 405 x 850mm

Code	Description
MH1103-0000	Trolley without trays
MH1103-0040	4 x 100mm trays
MH1103-0031	3 x 100mm trays    1 x 200mm tray
MH1103-0012	1 x 100mm tray    2 x 200mm trays
MH1103-0230	2 x 50mm trays    3 x 100mm trays



MH1103-0040 shown

### LOW LEVEL EURO SINGLE - 668 x 405 x 850mm

Code	Description
MH1106-0000	Trolley without trays
MH1106-0040	4 x 100mm trays
MH1106-0031	3 x 100mm trays    1 x 200mm tray
MH1106-0012	1 x 100mm trays    2 x 200mm trays
MH1106-0230	2 x 50mm trays    3 x 100mm trays



MH1106-0040 shown

### LOW LEVEL COMPACT DOUBLE - 911 x 405 x 850mm

Code	Description
MH1203-0000	Trolley without trays
MH1203-0080	8 x 100mm trays
MH1203-0062	6 x 100mm trays    2 x 200mm trays
MH1203-0024	2 x 100mm trays    4 x 200mm trays
MH1203-0460	4 x 50mm trays    6 x 100mm trays



MH1203-062 shown

### LOW LEVEL EURO DOUBLE - 911 x 600 x 850mm

Code	Description
MH1206-0000	Trolley without trays
MH1206-0080	8 x 100mm trays
MH1206-0062	6 x 100mm trays    2 x 200mm trays
MH1206-0024	2 x 100mm trays    4 x 200mm trays
MH1206-0460	4 x 50mm trays    6 x 100mm trays



MH1206-0460 shown





Ideal for everyday use on busy hospital wards



Fixed rails on both side provide balanced handling

- One piece fully welded frame
- Choice of tray configurations
- Trays supplied with dividers
- 100mm castors and 2 braked
- Four ring buffers
- Removable sides top and back
- Easy cleaning system
- 10 year guarantee

## Standard height with handles - 1050mm high

**OPTIONAL ACCESSORIES SEE PAGE 7**

**FOR PETO ONLINE ORDERING SYSTEM:**

Please add VO1ZZ to the start of the product code.  
E.G. VO1ZZMH3103-0000

**COMPACT TRAY 400 x 300mm CONTAINS:**

1 lengthway and widthway divider.  
2 label holders and 10 coloured inserts in 5 colours.

**EURO TRAY 600 x 400mm CONTAINS:**

2 lengthway and widthway dividers.  
2 label holders and 10 coloured inserts in 5 colours.

**FRAME COLOUR OPTIONS:** Lilac available from stock



Lilac



Pale green



White



Grey



Red



Lilac frame colour is available from stock for fast track delivery

### STANDARD HEIGHT WITH HANDLES COMPACT SINGLE 468 x 405 x 1050mm

Code	Description		
MH3103-0000	Trolley only		
MH3103-0050	5 x 100mm trays		
MH3103-0031	3 x 100mm trays	1 x 200mm tray	
MH3103-0003	3 x 200mm trays		
MH3103-0112	1 x 50mm tray	1 x 100mm tray	2 x 200mm trays
MH3103-0320	3 x 50mm trays	2 x 100mm trays	



MH3103-0031 shown

### STANDARD HEIGHT WITH HANDLES SINGLE EURO 668 x 405 x 1050mm

Code	Description		
MH3106-0000	Trolley only		
MH3106-0050	5 x 100mm trays		
MH3106-0031	3 x 100mm trays	1 x 200mm trays	
MH3106-0003	3 x 200mm trays		
MH3106-0112	1 x 50mm trays	1 x 100mm trays	2 x 200mm trays
MH3106-0320	3 x 50mm trays	2 x 100mm trays	



MH3106-0112 shown

### STANDARD HEIGHT WITH HANDLES COMPACT DOUBLE 911 x 405 x 1050mm

Code	Description		
MH3203-0000	Trolley only		
MH3203-0010	10 x 100mm trays		
MH3203-0062	6 x 100mm trays	2 x 200mm trays	
MH3203-0006	6 x 200mm trays		
MH3203-0224	2 x 50mm trays	2 x 100mm trays	4 x 200mm trays
MH3203-0640	6 x 50mm trays	4 x 100mm trays	



MH3203-0062 shown

### STANDARD HEIGHT WITH HANDLES EURO DOUBLE 911 x 600 x 1050mm

Code	Description		
MH3206-0000	Trolley only		
MH3206-0010	10 x 100mm trays		
MH3206-0062	6 x 100mm trays	2 x 200mm trays	
MH3206-0006	6 x 200mm trays		
MH3206-0224	2 x 50mm trays	2 x 100mm trays	4 x 200mm trays
MH3206-0640	6 x 50mm trays	4 x 100mm trays	



MH3206-0010 shown

## Standard height - 950mm high

One piece fully welded frame  
Shatterproof trays  
10 year guarantee

### OPTIONAL ACCESSORIES SEE PAGE 7

#### FOR PETO ONLINE ORDERING SYSTEM:

Please add VO1ZZ to the start of the product code.  
E.G. VO1ZZMH2103-0000

#### COMPACT TRAY 400 x 300mm CONTAINS:

1 lengthway and widthway divider.  
2 label holders and 10 coloured inserts in 5 colours.

#### EURO TRAY 600 x 400mm CONTAINS:

2 lengthway and widthway dividers.  
2 label holders and 10 coloured inserts in 5 colours.

#### FRAME COLOUR OPTIONS: Lilac available from stock



Lilac Pale green White Grey Red



Lilac frame colour is available from stock for fast track delivery

### STANDARD HEIGHT COMPACT SINGLE 468 x 405 x 950mm

Code	Description
MH2103-0000	Trolley only
MH2103-0050	5 x 100mm trays
MH2103-0031	3 x 100mm trays 1 x 200mm tray
MH2103-0003	3 x 200mm trays
MH2103-0112	1 x 50mm tray 1 x 100mm tray 2 x 200mm trays
MH2103-0320	3 x 50mm trays 2 x 100mm trays



MH2103-0031 Shown

### STANDARD HEIGHT SINGLE EURO 668 x 405 x 950mm

Code	Description
MH2106-0000	Trolley only
MH2106-0050	5 x 100mm trays
MH2106-0031	3 x 100mm trays 1 x 200mm tray
MH2106-0003	3 x 200mm trays
MH2106-0112	1 x 50mm trays 1 x 100mm tray 2 x 200mm trays
MH2106-0320	3 x 50mm trays 2 x 100mm trays



MH2106-0050 shown

### STANDARD HEIGHT COMPACT DOUBLE 911 x 405 x 950mm

Code	Description
MH2203-0000	Trolley only
MH2203-0010	10 x 100mm trays
MH2203-0062	6 x 100mm trays 2 x 200mm trays
MH2203-0006	6 x 200mm trays
MH2203-0224	2 x 50mm trays 2 x 100mm trays 4 x 200mm trays
MH2203-0640	6 x 50mm trays 4 x 100mm trays



MH2203-0062 shown

### STANDARD HEIGHT EURO DOUBLE 911 x 600 x 950mm

Code	Description
MH2206-0000	Trolley only
MH2206-0010	10 x 100mm trays
MH2206-0062	6 x 100mm trays 2 x 200mm trays
MH2206-0006	6 x 200mm trays
MH2206-0224	2 x 50mm trays 2 x 100mm trays 4 x 200mm trays
MH2206-0640	6 x 50mm trays 4 x 100mm trays



MH2206-0062 shown



## Tall trolleys - 1450mm high

High capacity mobile storage units  
Easy to manoeuvre  
Simple to clean  
10 year guarantee

**FRAME COLOUR OPTIONS:** Lilac available from stock



Lilac    Pale green    White    Grey    Red

**OPTIONAL ACCESSORIES SEE PAGE 7**

**FOR PETO ONLINE ORDERING SYSTEM:**

Please add VO1ZZ to the start of the product code.  
E.G. VO1ZZMH4106-0000

**COMPACT TRAY 400 x 300mm CONTAINS:**

1 lengthway and widthway divider.  
2 label holders and 10 coloured inserts in 5 colours.

**EURO TRAY 600 x 400mm CONTAINS:**

2 lengthway and widthway dividers.  
2 label holders and 10 coloured inserts in 5 colours.

### TALL TROLLEY SINGLE EURO 668 x 405 x 1450mm

Code	Description
MH4106-0000	Trolley only
MH4106-0080	8 x 100mm trays
MH4106-0033	3 x 100mm trays    3 x 200mm trays
MH4106-0014	1 x 100mm trays    4 x 200mm trays
MH4106-0142	1 x 50mm trays    4 x 100mm trays    2 x 200mm trays
MH4106-0052	5 x 100mm trays    2 x 200mm trays
MH4106-0005	5 x 200mm trays



MH4106-0080 shown

### TALL TROLLEY DOUBLE EURO 911 x 600 x 1450mm

Code	Description
MH4206-0000	Trolley only
MH4206-0916	16 x 100mm trays
MH4206-0066	6 x 100mm trays    6 x 200mm trays
MH4206-0028	2 x 100mm trays    8 x 200mm trays
MH4206-0284	2 x 50mm trays    8 x 100mm trays    4 x 200mm trays
MH4206-0014	10 x 100mm trays    4 x 200mm trays
MH4206-0001	10 x 200mm trays



MH4206-0284 shown

## Stainless steel dressing & Instrument trolleys

Hygienic seamless welded frame

304 stainless steel construction

Fixed or removable shelves

4 x 75mm anti-static castors



**ACCESSORIES:**

Customise the trolley design by suffixing the accessories codes.

E.G: MS4401BBC

Dressing trolley, protective buffers and braked castors.



**Code Description**

- B Protective buffers
- BC Braked castors
- SD Single drawer
- SDL Single drawer with lock
- DD Double drawer
- DDL Double drawer with lock
- GR Guide rail
- IV IV pole

**FOR PETO ONLINE ORDERING SYSTEM:**

Please add VO1ZZ to the start of the product code.

E.G. VO1ZZMS4401

**STAINLESS STEEL DRESSING TROLLEY 460 x 460 x 890mm**

Code	Description
MS4401	Removable shelves 3 x flanges up 1 x flange down
MS4402	Fixed shelves 3 x flanges up 1 x flange down
MS4403	Removable shelves Flanges down
MS4404	Fixed shelves Flanges down
MS4405	Removable shelves Flanges up



MS4401B shown with buffers

**STAINLESS STEEL DRESSING TROLLEY 610 x 460 x 890mm**

Code	Description
MS6401	Removable shelves 3 x flanges up 1 x flange down
MS6402	Fixed shelves 3 x flanges up 1 x flange down
MS6403	Removable shelves Flanges down
MS6404	Fixed shelves Flanges down
MS6405	Removable shelves Flanges up



MS6401 shown

**STAINLESS STEEL DRESSING TROLLEY 920 x 460 x 890mm**

Code	Description
MS9402	Fixed shelves 3 x flanges up 1 x flange down
MS9404	Fixed shelves Flanges down



MS9402 shown

**STAINLESS STEEL DRESSING TROLLEY 1220 x 460 x 890mm**

Code	Description
MS1402	Fixed shelves 3 x flanges up 1 x flange down
MS1404	Fixed shelves Flanges down



MS1402 shown

**STAINLESS STEEL DRESSING TROLLEY 1220 x 610 x 890mm**

Code	Description
MS1602	Fixed shelves 3 x flanges up 1 x flange down
MS1604	Fixed shelves Flanges down



MS1602 shown



## I.V. poles and drip stands

CE Marked stands

Class 1 devices

Conform to BS3619:1976

- Five leg design with low centre of gravity to assist stability
- Heavy duty base with plump rounded corners for protection of the users, patients & fabric of the building.
- Fully welded construction with 30mm diameter scratch resistant stainless steel upright.
- Polymer 50mm twin wheel castors shielded from fluid spills with optional breaks.
- Height adjustable stainless steel pole using a handwheel (1200 to 2000mm)
- Available with two/four seamlessly welded pigtail hooks
- Maximum hook load 4kg (per stand)
- Suitable for use with two infusion pumps/syringe drivers (max 4kg per stand).

### FOR PETO ONLINE ORDERING SYSTEM:

Please add VO1ZZ to the start of the product code.  
E.G. VO1ZZMS3302


### TWO HOOK I.V. POLE/DRIP STAND 200 - 1200mm

Code	Description
MS3302	Two hook stand
MS3302-WB	Two hook stand with weighted base pole

### FOUR HOOK I.V. POLE/DRIP STAND 200 - 1200mm

Code	Description
MS3304	Four hook stand
MS3304-WB	Four hook stand with weighted base pole

### ANGLE TILT TEST



### Ref Description

- A** Test 1 stability test to BS 3619:1976 with single infusion 1000ml dispenser 9° permissible angle achieved.
- B** Test 2 with Flo-gard pump at 0.8m height and infusion dispenser at 1.9m 11° permissible angle achieved.
- C** Test 3 with Flo-gard pump at 1.0m height and infusion dispenser at 1.9m 9° permissible angle achieved.



MS3302 shown



Two hook design



Weighted base pole



MS3304 shown



Four hook design



Weighted base pole



## Recharge trolley

Charges from one mains outlet  
 Single or double sized unit  
 Mobile unit with castors  
 Full CE safety certification


**FRAME COLOUR OPTIONS:** Lilac available from stock



### FOR PETO ONLINE ORDERING SYSTEM:

Please add VO1ZZ to the start of the product code.  
 E.G. VO1ZZMJ4450-1005

### FEATURES:



MJ4450-1005 shown

Ref	Description
	<p><b>A</b> One piece mild steel welded frame with a selection of colour options. No site assembly required.</p> <p><b>B</b> Fitted with 4 anti-static 100mm castors - Non marking, buffered wheels with locking brakes.</p> <p><b>C</b> UK mains sockets 230v 50 Hz            3 x double sockets on single unit            6 x double sockets on double unit</p> <p><b>D</b> RCD safety cut out</p> <p><b>E</b> Emergency stop button</p> <p><b>F</b> Medstor modular adjustable longspan shelves with 4 hook in support brackets.</p>

Single recharge trolley 1040 x 462 x 1435mm

Code	Description
MJ4450-1005	Single unit with 6 sockets

Double recharge trolley 2030 x 462 x 1435mm

Code	Description
MJ4450-2005	Double unit with 12 sockets



Lilac case notes trolley with file holder



Medical records transfer trolley

## Medical records and Case note trolleys

One piece robust construction  
Added enhanced functionalities

### FOR PETO ONLINE ORDERING SYSTEM:

Please add VO1ZZ to the start of the product code.  
E.G. VO1ZZMS3001

**FRAME COLOUR OPTIONS:** Lilac available from stock



Lilac Pale green White Grey Red



### FILING POCKETS:

Ref Code Description

Ref Code	Description
A MT4480-5002	Double filing pocket combined X-ray Foolscape 350mm wide x 460mm high Material Faux leather
B MT4480-5003	Filing pocket X-ray 350mm wide x 460mm high Material Faux leather
C MT4480-5004	Filing pocket Foolscape 350mm wide x 280mm high Material Faux leather

### CASE NOTE TROLLEYS LOCKABLE

- One piece welded frame, coloured coated mild steel.
- Enclosed suspension filing area and drop down laminated worktop.
- Integral sheet storage system.
- Lower shelf for larger note file.
- Fixed recessed base shelf and external storage pockets.
- 4 x 100mm castors.
- Optional Buffers.

Code	Description
MS3001	Medium case notes trolley 975 x 475 x 980mm
MS3100	Small case notes trolley 625 x 475 x 980mm



Rails fitted for filing pockets

MS3100 shown with MM4479FP file holder

### MEDICAL RECORDS TRANSFER

- One piece welded frame, coloured coated mild steel.
- Enclosed suspension filing area. Lockable
- Fixed recessed base shelf
- 4 x 100mm castors.
- Optional Buffers

Code	Description
MS3002	Medical records transfer trolley 840 x 550 x 980mm Box size 610 x 470 x 233mm



Secure locking for privacy

MS3002 shown



Light duty general purpose trolley



General purpose MS3104 in lilac

## General purpose

Stainless steel, mild steel,  
Polymer and utility trolleys

### FOR PETO ONLINE ORDERING SYSTEM:

Please add VO1ZZ to the start of the product code.  
E.G. VO1ZZMS3103HD

**FRAME COLOUR OPTIONS:** Lilac available from stock



Lilac



Pale green



White



Grey



Red

### GENERAL PURPOSE STAINLESS STEEL

- One piece frame construction.
- Integral push handles.
- Recessed shelves to contain spills
- 4 x 100mm castors..

Code	Description
MS3103-HD	Heavy duty Stainless steel trolley 1260 x 485 x 1050mm
MS3103-LD	Light duty Stainless steel trolley 1260 x 485 x 1050mm



MS3103HD shown

### LIGHT DUTY GENERAL PURPOSE THREE TIER

- Manufactured from tubular steel
- Choice of colour finishes
- Drop in moisture resistant shelves
- 4 x 100mm castors.
- Optional Buffers

Code	Description
MS3104	Medical records transfer trolley 1260 x 485 x 1050mm



MS3104 shown with optional buffers



## Linen and waste

75% saving against plastic laundry bags  
Boiling and tumble-dry qualities  
1000 wash guarantee

### FOR PETO ONLINE ORDERING SYSTEM:

Please add VO1ZZ to the start of the product code.  
E.G. VO1ZZMS12ASM

### BAG COLOUR OPTIONS: PLEASE SPECIFY AT POINT OF ORDER



Blue



Green



White



Grey



Red

Ensuring soiled and infected laundry items do not come into contact during transportation to the washing machine.

The bags automatically unwrap themselves and release the laundry into the machine, avoiding the spread of germs.

Code	Description	Size
MS12ASM	Linen bag 100% polyester	460 x 430 x 870mm

### LINEN SINGLE BAG HOLDER

- Sprung frame linen bag holders
- Use with fabric or polythene bags
- Easy clean panels
- Non marking castors

#### FOR PETO ONLINE ORDERING SYSTEM:

Please add VO1ZZ to the start of the product code.  
E.G. VO1ZZMS3206

#### SUITABLE FOR:

Fabric or polythene bags 1500mm circumference.  
Adjustable to suit other bag sizes. Minimum bag height 750mm

Code	Description	Size
MS3206	Single bag holder	460 x 400 x 825mm

#### FRAME COLOUR OPTIONS: Lilac available from stock



MS3206 shown with linen bag

### MILD STEEL LINEN BAG HOLDERS

- Steel linen bag trolley
- For use with 1500 circumference bags
- Easy clean panels
- Coloured mild steel

#### SUITABLE FOR:

Fabric or polythene bags 1500mm circumference.  
Adjustable to suit other bag sizes. Minimum bag height 750mm

Code	Description	Size
MS3201PC	Single bag holder	550 x 380 x 1060mm
MS3202PC	Double bag holder	880 x 420 x 910mm

#### FRAME COLOUR OPTIONS: Lilac available from stock



MS3201PC shown



MS3202PC shown with linen bag

### STAINLESS STEEL LINEN BAG HOLDERS

- Steel linen bag trolley
- For use with 1500 circumference bags
- Fitted with 125mm castors

Code	Description	Size
MS3201SS	Single bag holder	550 x 380 x 1060mm
MS3202SS	Double bag holder	880 x 420 x 910mm

#### FOR PETO ONLINE ORDERING SYSTEM:

Please add VO1ZZ to the start of the product code.  
E.G. VO1ZZMS3201SS



MS3201SS shown



MS3202SS shown

### LINEN SINGLE BAG HOLDER WITH LID

- Linen bag trolley with lid
- Pedal operated aluminium formed lid
- 2 fixed and 2 swivel castors
- Selection of lid colours

#### FRAME COLOUR OPTIONS: Lilac available from stock



Code	Description	Size
MS3200	Single bag holder	460 x 430 x 870mm

#### FOR PETO ONLINE ORDERING SYSTEM:

Please add VO1ZZ to the start of the product code.  
E.G. VO1ZZMS3200



MS3200 shown with and without bag

## Linen and waste

10 year guarantee  
Welded frame construction

### FOR PETO ONLINE ORDERING SYSTEM:

Please add VO1ZZ to the start of the product code.  
E.G. VO1ZZMS3205

**FRAME COLOUR OPTIONS:** Lilac available from stock



Lilac



Pale green



White



Grey



Red

### HAMPER STYLE LINEN TROLLEY

- Four secure corner fixings
- One piece welded frame for strength
- 100mm non marking castors

Code	Description	Size
MS3205	4 corner fixing - 200kg load*	820 x 910 x 995mm

\* Max working load recommended



MS3205 is fitted with 4 x 100mm castors

### COMBINATION LINEN TROLLEY

- Central panel for segregation
- Removable aluminium trays
- 10 year guarantee

Code	Description	Size
MS3208	Bag holder with 3 trays	815 x 483 x 1130mm



MS3208 trolley shown with linen bag

### MOBILE LINEN TROLLEY

- Steel framework with aluminium panels
- 2 adjustable shelves
- 4 buffered castors with 2 front locking
- Available with protective fabric front cover



Protective front cover



MS3207 shown

Code	Description	Size
MS3207	Storage trolley	1170 x 510 x 1030mm
MS3207-FC	Pull down protective fabric front cover for MS3207 Non Velcro fastening	

### FIRE RETARDANT SACK HOLDER

- Manufactured from high density polyethylene
- Tested to NHS technical memorandum 83 (Fire safety)
- All units fitted with bag cord
- Moulded-in graphics available additional cost
- Stainless steel push rod
- Lid does not touch the wall when opening

**LID COLOURS:** Please suffix code:

Suffix B



Black

Suffix O



Orange

Suffix Y



Yellow



Pedal operated



Silent closing lid option\*

Code	Description	Size
MT0036	20 litre fixed sack holder	620 x 240 x 405mm
MT0037	65 litre mobile sack holder	810 x 400 x 480mm



MT0036 shown



MT0037 shown



For use with wet-rooms and cleaning areas



Lilac cylinder trolleys are available via fast track delivery

## Cylinder holders

Compatible with a selection of cylinders  
Secure fittings

### FOR PETO ONLINE ORDERING SYSTEM:

Please add VO1ZZ to the start of the product code.  
E.G. VO1ZZMS3502

### CYLINDER TROLLEYS

- Holds cylinders securely in place
- Captive chain design
- Anti static wheels

Code	Description	Size
MS3502	A/F sized holder	305 x 265 x 820mm
MS3503	F or G type holder	560 x 420 x 1150mm

FRAME COLOUR OPTIONS: Lilac available from stock



MS3502 shown



MS3503 shown

### PORTABLE OXYGEN CARRIER

- Balanced carry handle
- Fits to most bed sides and medi-rails
- Wall bracket for portable cylinder holder

Code	Description	Size
MS3504	D & E sized cylinders	390 x 150 x 385mm
MS3501	Fits D cylinders	880 x 420 x 910mm
MS3501-WB	Wall bracket for MS3501	



MS3504 shown



MS3501 shown with wall bracket

### DOUBLE CYLINDER WALL BRACKET HOLDER

- One piece design
- Captive chain to secure cylinders
- Wall fixing pre-drilled holes
- Mild steel

Code	Description	Size
MS3500-01	Fits D or E cylinders	359 x 107 x 88mm
MS3500-02	Fits F cylinders	435 x 126 x 88mm
MS3500-03	Fits G cylinders	511 x 145 x 88mm
MS3500-04	Fits J cylinders	613 x 170 x 88mm



MS3500-01 shown



Fits variety of sizes

# MedStor<sup>®</sup>

The heart of your storage solutions



**Mailbox<sup>®</sup>**  
Stamford Products Limited

**UK:**

Bayley Street, Stalybridge, Cheshire, SK15 1QQ.

Tel: 0161 330 5577 Fax: 0161 330 5576

Website: [www.medstor.co.uk](http://www.medstor.co.uk) email: [health@medstor.co.uk](mailto:health@medstor.co.uk)

**Ireland:**

Unit 120, Grange Hill, Baldoye Industrial Estate, Dublin 13.

Tel: 01 839 3596 Fax: 01 839 3468

Website: [www.medstor.ie](http://www.medstor.ie) email: [sales@medstor.ie](mailto:sales@medstor.ie)